



GMP Power Tug

GMP Power Tug:

The EHS Solutions GMP Power Tug is designed to dock with various heavy mobile equipment, ensuring safe and easy transportation throughout process areas.

Max Weight Capacity:

- Standard weight class option: 6000lbs

Drive System (with utility requirements):

- Battery-Powered
 - 24-volt 65AH maintenance-free batteries
 - Integrated 110VAC plug-in automatic charger
- 24-Volt permanent magnet DC drive motor
- Electromagnetic disc brake
- On/Off key switch and E-stop

Material of Construction Options:

- 304 Stainless Steel Shrouding
- 316 Stainless Steel Shrouding

Surface Finish Grade (see Section 2c: Surface Finish Grade Options):

- GMP Standard (Random Orbital Polish; 120 micro-inch Ra)
- GMP + (Unidirectional Polish; 35 micro-inch Ra)

Power Tug Attachments

- Custom Per Application



POWER TUG DETAILS

