

## soloBAG

### Keep a tighter grip on costs with our glove bags

Glove bags are at the heart of what we do – so we've made sure our products have significant price AND performance advantages over our competition.



soloBAG is our tag for Glove Bags: soloBAGs can be used to enclose any process, maintenance or clinical compounding application where the operation calls for a short life cycle. These entry-level polyurethane/PVC flexible envelopes are supplied with a simple breather filter to permit contamination free envelope collapse at the end of the life cycle. For single use glove bags and connection sleeves Solo can offer fully custom designed LDPE glove bags to drive down costs even further.

### Superior design and better protection

Our glove bags offer an improved operator interface, reducing fatigue

Simple envelope mount systems

Integrated filter systems disposed of with the enclosure, reducing risk as there's no need to change filters at end of batch/campaign

Cross contamination risk eliminated by means of disposing the enclosure

Standard Operating Procedure (SOP) development to ensure safe operating cycle.

### Cost effective and easy to use

Reduced capital cost for a small batch/campaign

Cleaning validation not necessary when enclosure is disposed of

Quick to assemble, set-up and validate

Shorter project life cycle resulting in reduced project management costs

Little or no installation required.

### Wide range of designs and application history

As with all of our containment products, we can manufacture glove bags to a range of specifications including custom designs to your order. We can design, manufacture and deliver products quickly and efficiently, with spares available in days.

Our glove bags are designed for a wide range of processes such as **Sampling, Dispensing, Charging, Offloading and Filtering**. soloBAG products are currently used in numerous markets, **chemical production, pharmaceutical manufacture, contract research and manufacture and low grade nuclear waste transfer**.

We produce a range of accessories for soloBAG products including: connection sleeves (process, utility), couplings and 'zipper' type doors.

Safety and operator protection are top priorities so you can be sure soloBAG products meet international pharmaceutical industry product and toxicity bandings.

*[Talk to us about glove bags/flexible envelopes](#)*

To find out about prices and delivery, call **0800 292 2306** or **01625 344014** or [contact us by email](#).